


Freely configurable colon with insufflation and suction capabilities makes this colonoscopy model ideal for insertion technique training.



Main Features

- A water-soluble lubricant is used to replicate the actual feel of the colon's interior.
- The colon tube can be freely configured to suit training requirements.
- The anus has an open/close function that allows you to perform insufflation and suction.
- The colon tube can be removed and replaced with a new one.
- The colon tube is washable and can be kept clean for next usage.



Manufactured by: KYOTOKAGAKU CO.,LTD.

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.

OLYMPUS[®]

OLYMPUS MEDICAL SYSTEMS CORP.
Shinjuku Monolith, 2-3-1 Nishi-Shinjuku, Shinjuku-ku, Tokyo 163-0914, Japan

For a complete listing of
sales and distribution locations visit:
www.olympus.com